



G-EYE[®] Endoscope

YOUR GI ENDOSCOPE CAN
ACCOMPLISH SO MUCH MORE



- Increased diagnostic yield in colonoscopy
- Endoscope stabilization during intervention
- On-demand two-balloon endoscopy
- Unique Controlled Withdrawal™

ENDOSCOPY, THE SMART WAY

smart[○]
medical systems

The G-EYE® Endoscope comprises a standard endoscope, having a permanently integrated, reusable balloon at its bending section. This novel balloon-endoscope enables increased diagnostic yield during colonoscopy, endoscope stabilization during therapy, and extended operation range in the colon, small bowel and entire GI tract. Employing a unique Controlled Withdrawal™ technique with the balloon partially inflated, the G-EYE® endoscope provides austral folds straightening and enhanced detection capabilities. Upon need, the balloon can be readily inflated to anchor and stabilize the

endoscope during intervention. The G-EYE® endoscope, combined with the NaviAid™ AB or NaviAid™ ABC through-the-channel balloon-devices, provides an ultimate assembly for performing deep and swift two-balloon endoscopy, while avoiding limitations of other balloon techniques. With the G-EYE® system there is no need for over-tube usage, mounting and pre-procedure preparations, or extended learning process. The G-EYE® balloon and the balloon of the AB or ABC devices are controlled simultaneously by the NaviAid™ SPARK²C inflation system.

ADDITIONAL KEY ADVANTAGES:

- Permanently integrated, reusable balloon
- No over-tube, no preparations prior to the procedure
- Both colonoscopy and upper GI endoscopy
- Intuitive, easy to operate, short learning curve
- Supports additional NaviAid™ devices

G-EYE® – PRODUCT SPECIFICATIONS:

Suitable Endoscopes	Various standard brands and models
Balloon Diameter	Up to 60mm
Nominal Anchoring Pressure	70mbar
Controlled Withdrawal™ Pressure	3 intermediate pressure levels
SPARK ² C Inflation System – Power Supply	100 - 240VAC; 50-60Hz



ABOUT SMART MEDICAL

Your professional partner in a variety of endoscopy applications, SMART is a pioneer in the development and manufacture of world-leading medical devices. SMART's NaviAid™ product family enhances the capabilities and range of standard endoscopes in an intuitive, cost-effective manner. SMART is headquartered in Israel's hi-tech region, and operates in the United States via its wholly-owned subsidiary, SMART GI Inc. SMART is ISO 13485 certified, and its globally-distributed commercial products have gained FDA clearance (in the USA) and the CE Mark (in Europe).



ENDOSCOPY, THE SMART WAY

Headquarters: 5 HaNofar st, P.O Box 4123 Raanana, 4366404 Israel. Tel:+972-9-7444321, Fax:+972-9-7444543, info@smartmedsys.com • US Offices: SMART GI Inc., 540 Frontage Road, Suite 2265, Northfield, IL 60093, Tel: 1-847-920-5442, Fax: 1-800-230-5823 us@smartmedsys.com